

Generating new synthesis PET images through out deep learning model (Topic -11)

Nguyen Thi Van,

MSc candidate student

**Department of Artificial Intelligence Convergence
Chonnam National University**

Student ID: 217366

Video Demo

The screenshot shows a Visual Studio Code (VS Code) interface with the following details:

- File Explorer:** On the left, it lists several Python extensions:
 - @category:debuggers
 - Python (IntelliSense (Pylance))
 - Microsoft
 - Python Preview (67ms)
 - AREPL for ... (551K)
 - real-time python scra...
 - Almenon (Install)
 - Python Image (239K)
 - Numpy, Pillow, Open...
 - 윤대희 (Install)
 - Python C++ (165K)
 - Extension for starting...
 - BeniBenj (Install)
 - GeeCode Python (77K)
 - Linting, Debugging (...)
 - GeeCode Team (Install)
 - LiveCode for ... (68K)
 - Real-time python vari...
 - xirider (Install)
 - Python Resource (42K)
 - A resource monitor fo...
 - kaih2o (Install)
 - Python Quick (46K)
 - Quickly print out log ...
 - AhadCove (Install)
 - Python Debug (32K)
 - Prompt for host/port ...
 - Alexander Sa... (Install)
 - Quick-Python (9K)
 - Quickly handle print o...
 - WeidaWang (Install)
 - Python Quick (3K)
 - Wrap variables with s...
 - Brian Crant (Install)
 - NI Python Data (8K)
 - Create and export NI ...
 - NI (Install)
- Editor:** The main editor area displays a Python script named `main_liver.py`. The code uses Streamlit to create a user interface for liver segmentation. It includes imports for Streamlit (`st`), pandas (`pd`), and other libraries. A function `liver_show()` is defined to handle patient information and X-ray processing.
- Terminal:** At the bottom, the terminal window shows the output of running the Streamlit application:

```
The default interactive shell is now zsh.  
To update your account to use zsh, please run `chsh -s /bin/zsh`.  
For more details, please visit https://support.apple.com/kb/HT208050.  
(base) imieis-MacBook-Pro:streamlit_liver imei$
```
- Bottom Status Bar:** The status bar at the bottom indicates the current line and column (Ln 138, Col 75), the number of spaces (Spaces: 4), the encoding (UTF-8), and the file type (LF Python).


Summarize for this project



1. Motivation and target of project




Transferring via creating new synthesis images between PET and CT images



2.Purpose



Supporting when combined modality method for diagnosis



Synthesis CT



Synthesis PET

3. Datasets





Private datasets
(from hospital with 30 patients)

Public datasets (from
MICCAI conference
with 900 patients)

Training from 60% public datasets;
Testing from 40% public datasets
and all private source

4.Pre-processing datasets



- 1.Hounsfield Units
2. Contrast Limited Adaptive Histogram Equalization (CLAHE)

4. Training model




Condition
GANs model

5. Metrics for valuate results





SSIM & MSE
Loss-function


6. GUI



4. Training model









5. Analyzing datasets



CT whole body



PET whole body




PET whole body showing via SUV values




CT axial



PET axial




PET axial via SUV values



CT sagittal



PET sagittal




PET sagittal via SUV values

After training


CT image



CT



PET U-Net GAN




PET Res-Net GAN




PET Ground Truth

PET image



PET



CT U-Net GAN



CT Res-Net GAN



CT Ground Truth


Axial




PET U-Net GAN



PET Res-Net GAN



PET Ground Truth




PET U-Net GAN



PET Res-Net GAN



PET Ground Truth



PET Ground Truth


PET U-Net GAN

PET Res-Net GAN

Axial

Coronal

Sagittal



PET Ground Truth

PET U-Net GAN

PET Res-Net GAN

Axial


Coronal

Sagittal

6. Results



CT



PET U-Net GAN




PET Res-Net GAN




PET Ground Truth


CT images synthesis with **whole body** images




PET



CT U-Net GAN




CT Res-Net GAN




CT Ground Truth

CT images synthesis with each slice for each position


1. Axial



2. Sagittal



3. Coronal



PET images synthesis with whole body images



PET ground truth



CT



PET GAN



PET ground truth




CT




PET GAN

Loss function with U-Net backbone




Loss function with Res-Net backbone




SSIM & MSE

| | Data augmentation | | Data non-augmentation | |
|-----|-------------------|---------|-----------------------|---------|
| | U-Net | Res-Net | U-Net | Res-Net |
| | SSIM | 0.65 | 0.72 | 0.52 |
| MSE | 0.16 | 0.09 | 0.25 | 0.17 |

Loss function with U-Net backbone



Loss function with Res-Net backbone




SSIM & MSE

| | Data augmentation | | Data non-augmentation | |
|------|-------------------|---------|-----------------------|---------|
| | U-Net | Res-Net | U-Net | Res-Net |
| SSIM | 0.65 | 0.72 | 0.52 | 0.31 |
| MSE | 0.16 | 0.09 | 0.25 | 0.17 |


7. Evaluation

7.1 Loss function

Loss function with U-Net backbone



Loss function with Res-Net backbone



7.2 SSIM & MSE

The Structural Similarity Index (SSIM) metric extracts 3 key *features* from an image:




- Luminance
- Contrast
- Structure

The comparison the similarity between the two images is performed on the basis of these 3 features.


$$\text{SSIM}(\mathbf{x}, \mathbf{y}) = \frac{(2\mu_x\mu_y + C_1)(2\sigma_{xy} + C_2)}{(\mu_x^2 + \mu_y^2 + C_1)(\sigma_x^2 + \sigma_y^2 + C_2)}.$$

| | Data augmentation | | Data non-augmentation | |
|------|-------------------|---------|-----------------------|---------|
| | U-Net | Res-Net | U-Net | Res-Net |
| SSIM | 0.65 | 0.72 | 0.52 | 0.31 |
| MSE | 0.16 | 0.09 | 0.25 | 0.17 |

8. Interface to show



Showing the information of patients



- 1. This function allow to choose PET-ACT, PET-FDG or CT images**
- 2. Adding window size to adjust Hounsfield Units**

Choose one of the following

PET/CT images

Choose one of the following dataset

data/HJ19002 CT,PT

Choose one of the following

ACT

Center:

50 - +

Window size:

90 - +

Niter


5

Kappa

40

CT: Hounsfield windowing

0 503



Supplementation

Adding data public (MICCAI)

Pages / Wiki / Collections

...


A whole-body FDG-PET/CT dataset with manually annotated tumor lesions (FDG-PET-CT-Lesions)

Created by Brittney Camp, last modified by Tracy Nolan on Oct 04, 2022

Summary

Purpose: To provide an annotated data set of oncologic PET/CT studies for the development and training of machine learning methods and to help address the limited availability of publicly available high-quality training data for PET/CT image analysis projects. This data can also be used for machine learning challenges, which is exemplified in the autoPET MICCAI 2022 competition: <https://autopet.grand-challenge.org/>.

Data: The anonymized publication of data was approved by the local ethics committee and data protection officer. 501 consecutive whole body FDG-PET/CT data sets of patients with malignant lymphoma, melanoma and non small cell lung cancer (NSCLC) as well as 513 data sets without PET-positive malignant lesions (negative controls) examined between 2014 and 2018 at the University Hospital Tübingen were included. All examinations were acquired on a single, state-of-the-art PET/CT scanner (Siemens Biograph mCT). The imaging protocol consists of a diagnostic CT scan (mainly from skull base to mid-thigh level) with intravenous contrast enhancement in most cases, except for patients with contraindications. The following CT parameters were used: reference dose of 200 mAs, tube voltage of 120 kV, iterative reconstruction with a slice thickness of 2 - 3 mm. In addition, a whole-body FDG-PET scan was acquired 60 minutes after I.V. injection of 300-350 MBq 18F-FDG. PET data were reconstructed using an ordered-subset expectation maximization (OSEM) algorithm with 21 subsets and 2 iterations and a gaussian kernel of 2 mm and a matrix size of 400 x 400.



⌚ Training and test cohort

Training cases: 1,014 studies (900 patients)

Test cases (final evaluation): 200 studies

Test cases (preliminary evaluation): 5 studies

A case (training or test case) consists of one 3D whole body FDG-PET volume, one corresponding 3D whole body CT volume and one 3D binary mask of manually segmented tumor lesions on FDG-PET of the size of the PET volume. CT and PET were acquired simultaneously on a single PET/CT scanner in one session; thus PET and CT are anatomically aligned up to minor shifts due to physiological motion.

Training set

Training data consists of 1,014 studies acquired at the University Hospital Tübingen and is made publicly available on [TCIA](#) (as DICOM, NiFTI and HDF5 files).

Data is available on TCIA: DOI [10.7937/gkr-xv29](https://doi.org/10.7937/gkr-xv29)

Data pre-processing and structure


In a pre-processing step, the TCIA DICOM files are resampled (CT to PET imaging resolution, i.e. same matrix size) and normalized (PET converted to standardized uptake values; SUV).

For the challenge, the pre-processed data will be provided in NiFTI format. PET data is standardized by converting image units from activity counts to standardized uptake values (SUV). We recommend to use the resampled CT (`CTres.nii.gz`) and the PET in SUV (`SUV.nii.gz`). The mask (`SEG.nii.gz`) is binary with 1 indicating the lesion. The training and test database have the following structure.

NIFTI


```
|--- Patient 1
|   |--- Study 1
|   |   |--- SUV.nii.gz    (PET image in SUV)
|   |   |--- CTres.nii.gz (CT image resampled to PET)
|   |   |--- CT.nii.gz    (Original CT image)
|   |   |--- SEG.nii.gz    (Manual annotations of tumor lesions)
|   |   |--- PET.nii.gz    (Original PET image as activity counts)
|   |--- Study 2
|   |   |--- ...
|--- Patient 2
|   |--- ...
```

Data Augmentation




Data augmentation

Rotating





Cropping center





Up-creasing this volume datasets to 3-5 times. This process avoid bias in training stage

GUI tool document


Bar width to adjust Standard uptake values (SUV)



GUI tool document



The same with GUI tool for Bone. The interface when show liver segmentation also similar



GUI tool document

Choose one of the following

PET/CT images

Choose one of the following dataset

data/HJ19002 CT,PT

Choose one of the following

ACT

Center:

50

Window size:

90

Niter

1 5

Kappa

10 40

CT: Hounsfield windowing

0 503

Choose one of the following

PET/CT images

Choose one of the following dataset

data/HJ19002 CT,PT

Choose one of the following

ACT

Center:

50

Window size:


90

Niter

1 5

Kappa

10 40



Choose one of the following

PET/CT images

Choose one of the following dataset

data/HJ19002 CT,PT

Choose one of the following

ACT

Center:

50

Window size:


90

Niter

1 5

Kappa

10 40



Choose one of the following

PET/CT images

Choose one of the following dataset

data/HJ19002 CT,PT

Choose one of the following

ACT

Center:

50

Window size:


90

Niter

1 5

Kappa

10 40



GUI tool document

Choose one of the following

Segmentation mask ▾

Choose one of the following dataset

data/HJ19002 CT,PT ▾

Choose one of the following

ACT| ▾

Center:

40 - +

Window size:


70 - +

Niter

1

Kappa

50





File and environment to set up for GUI

Requirement

(1)


Detail package code



Library (Pip install these library in Python)

```
streamlit~=1.5.0
pydicom~=2.2.2
numpy~=1.22.1
ipywidgets~=7.6.5
Pillow~=9.0.0
ipython~=8.0.1
matplotlib~=3.5.1
ipyvolume~=0.5.2
sklearn~=0.0
scipy~=1.7.3
scikit-image~=0.19.1
plotly~=5.5.0
SimpleITK~=2.1.1
scikit-learn~=1.0.2
pandas
```

(2) Python version ≥ 3 .




Python 3.8.0

Release Date: Oct. 14, 2019

This is the stable release of Python 3.8.0.

Note: The release you're looking at is Python 3.8.0, an outdated release. Python 3.9 is now the latest feature release series of Python 3. Get the latest release of 3.9.x here.

(3) Anaconda



Link download:

<https://www.python.org/downloads/release/python-380/>

<https://www.anaconda.com/products/individual>

File and environment to set up for GUI

To detail for step to install and setup Python, Streamlit and Anaconda

Run package code on terminal by command line

1. Set up and install anaconda detail in this link: <https://docs.anaconda.com/anaconda/install/windows/>

2. Set up python detail in this link: <https://realpython.com/installing-python/>

3. Create enviroment "Streamlit": <https://docs.streamlit.io/library/get-started/installation>

3.1 Navigate to “streamlit_liver” folder and type command line on terminal as below

`cd streamlit_liver`

3.2. Create a new Pipenv environment in that folder and activate that environment by command line:

`pipenv shell`

3.3. Install Streamlit in your environment: `pip install streamlit`

3.4. Test that the installation worked: `streamlit hello`

The end!